

ASTOCOVER

Passive warming insulation

The ASTOCOVER can be used in pre-operative, peri-operative and post-operative situations, in intensive care or on the ward, as prophylaxis against the patient cooling down and as a treatment for hypothermia. With ASTOCOVER it is prevented that warmth is radiated from the patient into the area and it comes thus to no unnecessary calorific loss of the patient. Hence the patient is protected against cooling down caused by operation theatre ventilation and draught.

An ideal solution to every application can be found in the above-mentioned scenarios by selecting applied parts which can be combined.

In combination with our warming systems such as ASTOPAD SYSTEM (Patient Warmer), ASTOFLO and ASTOTHERM PLUS (Blood and Infusion Warmers) a complete conception is realized to prevent the patient from the negative side effects of an operation such as hypothermia and decubitus.



ASTOCOVER PAS100 Upper body cover
for passive warming insulation, 100 X 100 cm

Order no.: PAS 100

Cover can be washed and disinfected
for ASTOCOVER PAS100, 100 X 100 cm

Order no.: COV 40100



ASTOCOVER PAS150 Body cover
for passive warming insulation, 50 X 150 cm

Order no.: PAS 150

Cover can be washed and disinfected
for ASTOCOVER PAS150, 50 X 150 cm

Order no.: COV 40150



ASTOCOVER PAS151 Body cover
for passive warming insulation, 80 X 150 cm

Order no.: PAS 151

Cover can be washed and disinfected
for ASTOCOVER PAS151, 80 X 150 cm

Order no.: COV 40151



ASTOCOVER PAS155 Arm - Chest cover
for passive warming insulation, 50 X 150 cm

Order no.: PAS 155

Cover can be washed and disinfected
for ASTOCOVER PAS155, 50 X 150 cm

Order no.: COV 40155



ASTOCOVER PAS180 Body cover
for passive warming insulation, 80 X 180 cm

Order no.: PAS 180

Cover can be washed and disinfected
for ASTOCOVER PAS180, 80 X 180 cm

Order no.: COV 40180



ASTOCOVER Patient Cap
for passive warming insulation

Size **34** Order.-No.: PAS 034 Size **52** Order.-No.: PAS 052

Size **38** Order.-No.: PAS 038 Size **56** Order.-No.: PAS 056 Size **58** Order.-No.: PAS 058

ASTOCOVER PAS161 Leg Bag
for passive warming insulation

Order no.: PAS 161



Material: The cover of ASTOCOVER consists of a PU coated PES fabric.

Characteristics:

- breathe-active
- virus and bacteria impermeable
- blood, secretion and urine impermeable
- PVC free
- flexible and no formation of wrinkles
- Ökotex 100, II Certificate

Cleaning: For cleaning and disinfection usual, alcohol-based cleaning agents and disinfectants may be used. The washable and disinfectable covers for ASTOCOVER (COV40xxx) additionally may be washed at 95°C.



STIHLER ELECTRONIC

Stihler Electronic GmbH
Julius-Hoelder-Strasse 36 • 70597 Stuttgart • Germany
Telefon +49 (0) 711-720670 • Telefax +49 (0) 711-7206757
www.stihlerelectronic.de • E-Mail: info@stihlerelectronic.de